




Infrared Non-Contact Temple Thermometer

Model: HC700



- Telephoto temperature measurement distance within 10cm (3.94 inch)
- 1 second measurement
- Temple temperature
- Object temperature
- Auditory fever alarm
- 9 memories
- °C / °F switchable
- Backlight
- Auto shut-off
- Low battery indicator
- Self-diagnosis for malfunction



Model	Qty per carton	Carton volume
HC700	30 pcs	0.067 cbm/ctn
QTY per 20'	QTY per 40'	QTY per 40'HQ
12510 pcs	25050 pcs	29100 pcs

Storage and Transportation Condition
 -25°C ~ 55°C (-13°F ~ 131°F)
 ≤ 85% RH, 700~1060 hPa

Patent Nos :
 ZL 2014 3 0054262.7

Clinically Accurate

ASTM E1965-98:2009

This Temple thermometer has been clinically tested in several big teaching hospitals based on the protocol of ASTM E1965-98:2009 regulatory standard, covering sufficient febrile and normal body temperature subjects with satisfied clinical repeatability and measuring accuracy comparing to the concurrent oral temperature measurement reading.